

POCKET VISION SCREENER

95
874

distance
20
800
equivalent

2 8 4 3

point	jaeger	$\frac{20}{400}$
26	16	$\frac{20}{200}$
14	10	$\frac{20}{100}$
10	7	$\frac{20}{70}$
8	5	$\frac{20}{50}$
6	3	$\frac{20}{40}$
5	2	$\frac{20}{30}$
4	1	$\frac{20}{25}$
3	1+	$\frac{20}{20}$

6 3 8 E M E

8 7 4 5 E M W

6 3 9 2 5 M E E

4 2 8 3 6 5 W E M

3 7 4 2 5 8 E W E

9 3 7 8 2 6 W M E

4 2 8 7 3 9 E W E

Card is held in gook light 14 inches from eye. Record vision for each eye separately with and without glasses. Presbyopic patients should read thru bifocal segment. Check myopes with glasses only

PUPIL GAUGE (mm.)

